

Venalink MDS® Formats



**Monitored Dosage System:
The best service for your patients**





Venalink MDS®

Venalink MDS® is a reconditioned blister packaging, where the pharmacist organises the pills and/or capsules that a given patient takes, in order to indicate the day and time in which they should be administered, according to the prescribed dosage regimen, while providing the information necessary to do so correctly.

Advantages for pharmacies

- ✓ The best professional service
- ✓ Enhances the role of the pharmacist
- ✓ Improves disease control

Advantages for the patient

- ✓ Medication prepared by a pharmacist
- ✓ Error-free treatments with no excess stock
- ✓ Quality and safety guarantee



VENALINK MDS® Standard

Practical and handy design.



- ✓ **Compatible with semi-automation.**
- ✓ **Available in clear and amber blister packs (to decrease the risk of photo-instability).**
- ✓ **Available in Spanish, Portuguese, English, French, German and Catalan. Boxes containing 100, 250 and 500 units.**

VENALINK MDS® Landscape

Visually separated doses. Higher cavity capacity than that of standard formats. Ideal for patients who prefer this dosage guidance.

- ✓ **Available in Spanish. Boxes containing 100 and 250 units.**

NEW



VENALINK MDS® Monthly

Same format as the standard one, but for monthly preparations. Each MDS represents one daily dose. 4 variants: breakfast, lunch, dinner and free choice.

- ✓ **Available in Spanish and English. Boxes containing 500 units.**





VENALINK MDS® Qube

Visual and comfortable design. Large cavity capacity, which can easily accommodate multiple drugs. Large inner space to provide information. Ideal for patients with many drugs to recondition.

- ✓ Compatible with automation and semi-automation.
- ✓ Available in Spanish, Portuguese, Italian, English, French, German and Catalan. Boxes containing 100 and 250 units.



VENALINK MDS® Mini-Qube

Similar to Qube, but smaller.

- ✓ Available in Spanish. Boxes containing 100 and 250 units.



VENALINK MDS® 5 Dose

Ideal for patients with 5 daily doses or to check medication dosage, 'if needed'. Large inner space to provide information.

- ✓ Available in Spanish, Portuguese and English. Boxes containing 100 and 250 units.

VENALINK MDS®

FlexRX Light

The safe and discreet way to take your medication inside and outside your home. Ideal for polymedicated and active patients.


The large-capacity cavities are die-cut, hence allowing for easy separation by dose or day.

The cavity's rounded corners absorb energy better, if the medication is dropped, hence decreasing the risk of 'jumping'.



It does not require extra adhesive labels, as the information addressed to the patient is printed on the MDS label itself (which is compatible with any laser printer).

All of the cavities are identified with information regarding the medication they contain.

THIS FORMAT allows you to individually remove and label each dose.



- ✔ Requires software use.
- ✔ Compatible with automation and semi-automation.
- ✔ Available in Spanish. Boxes containing 100, 250 and 500 units.





We help you
help your
patients

More information:
+34 932 059 435 · venalink.es/en